



## Scholten Novatome™ Endomyocardial Biopsy Forceps

A biopsy design proven by over 45 years of successful use.



Sharp and steerable advantages in a disposable format.

Scholten Novatome™ is specially designed to acquire maximum specimen size for more accurate pathology results. The two-ringed handle easily adapt to any thumb position. The handle provides a solid instrument feel.

With a reputation for its ease of positioning and heart rhythm feedback, Scholten Bioptomes have earned high praise. The Novatome™ single use, disposable unit eliminates the risk of contamination. Scholten model is the best choice for safety and precision. Clearly, with a pedigree of more than 45 years of clinical experience, the disposable model Novatome™ will continue to set the standard worldwide.

Features:

- Perfect balance of flexibility and stiffness offering precise navigation and positioning
- Formable distal curve
- Familiar locking finger ring actuation handle
- Super precision machining
- Provides precise cutting edge alignment for every cut
- Single action cutting jaw (single moving jaw) provides superior cutting force
- Unique hollow ground cutting edge for superior sharpness



## Ordering Information

Ref no.	French size	Length	Curved tip	Jaw cup dia.	Qty/box
SU104-50	6	50 cm	Yes	2.0 mm	1
SU104-100	6	100 cm	No	2.0 mm	1
SU103-50	7	50 cm	Yes	2.3 mm	1
SU103-100	7	100 cm	No	2.3 mm	1
SU103-100 PC	7	100 cm	Yes	2.3 mm	1
SU102-50	8	50 cm	Yes	2.7 mm	1
SU101-50	9	50 cm	Yes	3.0 mm	1

MM 02-6000-05 ENG Rev. 02

### Contact Details, Nordics

customer@mermaidmedical.com  
mermaidmedical.com

### Contact Details, NL

info@mermaidmedical.nl  
mermaidmedical.com

### Contact Details, UK

customer@mermaidmedical.co.uk ·  
mermaidmedical.com