



Ultra-Pro 3™ Disposable Replacement Kits

The Ultra-Pro 3™ offers all the advantages of needle guidance while maintaining the “freehand” technique experience.

With zero needle resistance it provides smooth movement and navigation during puncture procedures.

The guide is designed for use with Ultra-Pro II™ and Ultra-Pro transducer brackets.

Features:

- Precision for reduced technique variability to rapidly deliver accuracy and stability for puncture procedures.
- Accurate positioning during puncture provides confidence for reaching the target.
- Two disposable inserts:
 - a large gauge insert offering 11-15 G.
 - a small gauge insert offering 16-22 G.



Quick-release feature for easy detachment of the needle from the transducer



Ultra-Pro-Series

Ordering Information

Ultra-Pro 3™ Disposable Replacement Kits

Ref no.	Descriptions	Qty/box
610-901	Sterile Ultra-Pro 3 with 14 x 91.5cm telescopically-folded CIV-Flex™ cover	24
610-919	Sterile Ultra-Pro 3 with 17.8 x 91.5cm telescopically-folded CIV-Flex™ cover	24

MM 02-1000-05 ENG Rev. 02.

Contact Details, Nordics

customer@mermaidmedical.com
mermaidmedical.com

Contact Details, NL

info@mermaidmedical.nl
mermaidmedical.com

Contact Details, UK

customer@mermaidmedical.co.uk ·
mermaidmedical.com