



Soft Tissue Biopsy Catalogue



SOFT TISSUE BIOPSY

Mermaid Medical offers a range of soft tissue biopsy needles designed to obtain high quality tissue samples with safety and efficiency in mind.

Our series of products include an automatic and semi-automatic option to support individual needs. In addition, a large range of gauge sizes, needle lengths, and throw lengths are available to match procedural requirements.

Mermaid Medical's range of soft tissue biopsy products includes diagnostic solutions for listed applications:



BREAST



KIDNEY



LIVER



LUNG



LYMPH



MUSCLE



PROSTATE



SPLEEN



THYROID

Mermaid Medical Group is a Nordic, privately owned company established in 2007 and headquartered in Copenhagen, Denmark. We develop, manufacture, and distribute medical devices to hospitals and end users across Europe, the U.S., and Asia. We primarily work within solutions to diseases in the vascular system as well as other devices used in interventional radiology. We strive to be the preferred partner for manufacturers as well as hospitals and healthcare professionals.

TABLE OF CONTENT



Automatic Disposable Biopsy Needle

- M-Biopsy® Automatic Disposable Biopsy Needle 4



Semi-automatic Disposable Biopsy Needle

- M-Biopsy® Semi-automatic Disposable Biopsy Needle 6






M·Biopsy® Automatic Disposable Biopsy Needle

The M·Biopsy® fully automatic disposable biopsy instrument is a user-friendly biopsy device designed to provide a high-quality biopsy sample.



It has been designed with focus on ergonomics to secure a comfortable grip, independent of device orientation. The left and right loading sliders secures a high-quality sample without compromising device fire power.

The M·Biopsy® coaxial needle comes with a sharp trocar tip and a transparent needle hub. The needle features an additional blunt stylet to minimize damage to tissue, vasculature, and organs.



Features:

- Ergonomic and lightweight one-handed design to secure comfortable grip, increase control, and minimize needle movement during imaging.
- 2-stage left and right loading button provides reduction of force needed to cock device.
- Comes with 11mm or 22mm throw length to meet a variety of clinical needs.
- Available as a single device or as a set with a coaxial needle for improved procurement and inventory control.






M·Biopsy® Semi-automatic Disposable Biopsy Needle

The M·Biopsy® semi-automatic biopsy instrument is a lightweight and adjustable device with the flexibility of a semi-automatic mode of action. The full loop scissor design provides a solid grip, leading to easy, single stage cocking and excellent maneuverability during the procedure. Color-coding on the inside of the handle loops makes identification of accurate gauge size easy.

In addition, its lightweight design combined with a wide range of gauge sizes and needle lengths also makes it ideal for CT-guided biopsies.

The M·Biopsy® coaxial needle comes with a sharp trocar tip and a transparent needle hub. The needle features an additional blunt stylet to minimize damage to tissue, vasculature, and organs.



Features:

- Adjustable throw length (10mm and 20mm) to accommodate individual procedural needs.
- Available as a single device or as a set with a coaxial needle for improved procurement and inventory control.
- Lightweight solid one-piece design makes it ideal for CT-guided biopsies.




Ordering and Specifications

M·Biopsy® Automatic Disposable Biopsy Needle

Ref. no.	Gauge	Length	Throw Length	Box	Ref. no.	Gauge	Length	Throw Length	Box
N60112100-CE	12	10cm	11mm	5	N70112100-CE	12	10cm	22mm	5
N60112130-CE	12	13cm	11mm	5	N70112130-CE	12	13cm	22mm	5
N60112160-CE	12	16cm	11mm	5	N70112160-CE	12	16cm	22mm	5
N60112200-CE	12	20cm	11mm	5	N70112200-CE	12	20cm	22mm	5
N60112250-CE	12	25cm	11mm	5	N70112250-CE	12	25cm	22mm	5
N60114100-CE	14	10cm	11mm	5	N70114100-CE	14	10cm	22mm	5
N60114130-CE	14	13cm	11mm	5	N70114130-CE	14	13cm	22mm	5
N60114160-CE	14	16cm	11mm	5	N70114160-CE	14	16cm	22mm	5
N60114200-CE	14	20cm	11mm	5	N70114200-CE	14	20cm	22mm	5
N60114250-CE	14	25cm	11mm	5	N70114250-CE	14	25cm	22mm	5
N60116100-CE	16	10cm	11mm	5	N70116100-CE	16	10cm	22mm	5
N60116130-CE	16	13cm	11mm	5	N70116130-CE	16	13cm	22mm	5
N60116160-CE	16	16cm	11mm	5	N70116160-CE	16	16cm	22mm	5
N60116200-CE	16	20cm	11mm	5	N70116200-CE	16	20cm	22mm	5
N60116250-CE	16	25cm	11mm	5	N70116250-CE	16	25cm	22mm	5
N60118100-CE	18	10cm	11mm	5	N70118100-CE	18	10cm	22mm	5
N60118130-CE	18	13cm	11mm	5	N70118130-CE	18	13cm	22mm	5
N60118160-CE	18	16cm	11mm	5	N70118160-CE	18	16cm	22mm	5
N60118200-CE	18	20cm	11mm	5	N70118200-CE	18	20cm	22mm	5
N60118250-CE	18	25cm	11mm	5	N70118250-CE	18	25cm	22mm	5
N60120100-CE	20	10cm	11mm	5	N70120100-CE	20	10cm	22mm	5
N60120130-CE	20	13cm	11mm	5	N70120130-CE	20	13cm	22mm	5
N60120160-CE	20	16cm	11mm	5	N70120160-CE	20	16cm	22mm	5
N60120200-CE	20	20cm	11mm	5	N70120200-CE	20	20cm	22mm	5
N60120250-CE	20	25cm	11mm	5	N70120250-CE	20	25cm	22mm	5

M·Biopsy® Automatic Disposable Biopsy Needle Set

Ref. no.	Gauge	Length	Throw Length	Matching Coaxial Needle			Box
				Gauge	Length	With Blunt Needle	
NS6011410-CE	14	10cm	11mm	13	5,9cm	Yes	5
NS6011413-CE	14	13cm	11mm	13	8,9cm	Yes	5
NS6011416-CE	14	16cm	11mm	13	11,9cm	Yes	5
NS6011610-CE	16	10cm	11mm	15	5,9cm	Yes	5
NS6011613-CE	16	13cm	11mm	15	8,9cm	Yes	5
NS6011616-CE	16	16cm	11mm	15	11,9cm	Yes	5
NS6011618-CE	16	18cm	11mm	15	13,9cm	Yes	5
NS6011810-CE	18	10cm	11mm	17	5,9cm	Yes	5
NS6011813-CE	18	13cm	11mm	17	8,9cm	Yes	5
NS6011816-CE	18	16cm	11mm	17	11,9cm	Yes	5
NS6011820-CE	18	20cm	11mm	17	15,9cm	Yes	5
NS6012010-CE	20	10cm	11mm	19	5,9cm	Yes	5
NS6012013-CE	20	13cm	11mm	19	8,9cm	Yes	5
NS6012016-CE	20	16cm	11mm	19	11,9cm	Yes	5
NS6012020-CE	20	20cm	11mm	19	15,9cm	Yes	5
NS7011410-CE	14	10cm	22mm	13	4,8cm	Yes	5
NS7011413-CE	14	13cm	22mm	13	7,8cm	Yes	5
NS7011416-CE	14	16cm	22mm	13	10,8cm	Yes	5
NS7011610-CE	16	10cm	22mm	15	4,8cm	Yes	5
NS7011613-CE	16	13cm	22mm	15	7,8cm	Yes	5
NS7011616-CE	16	16cm	22mm	15	10,8cm	Yes	5
NS7011810-CE	18	10cm	22mm	17	4,8cm	Yes	5
NS7011813-CE	18	13cm	22mm	17	7,8cm	Yes	5
NS7011816-CE	18	16cm	22mm	17	10,8cm	Yes	5
NS7011820-CE	18	20cm	22mm	17	14,8cm	Yes	5
NS7012010-CE	20	10cm	22mm	19	4,8cm	Yes	5
NS7012013-CE	20	13cm	22mm	19	7,8cm	Yes	5
NS7012016-CE	20	16cm	22mm	19	10,8cm	Yes	5
NS7012020-CE	20	20cm	22mm	19	14,8cm	Yes	5
NS7012025-CE	20	25cm	22mm	19	17,8cm	Yes	5




Ordering and Specifications

M·Biopsy® Semi-automatic Disposable Biopsy Needle

Ref. no.	Gauge	Length	Throw Length	Box
N50114060-CE	14	6cm	10/20mm	5
N50114100-CE	14	10cm	10/20mm	5
N50114130-CE	14	13cm	10/20mm	5
N50114160-CE	14	16cm	10/20mm	5
N50116060-CE	16	6cm	10/20mm	5
N50116100-CE	16	10cm	10/20mm	5
N50116130-CE	16	13cm	10/20mm	5
N50116160-CE	16	16cm	10/20mm	5
N50116200-CE	16	20cm	10/20mm	5
N50118060-CE	18	6cm	10/20mm	5
N50118100-CE	18	10cm	10/20mm	5
N50118130-CE	18	13cm	10/20mm	5
N50118160-CE	18	16cm	10/20mm	5
N50118200-CE	18	20cm	10/20mm	5
N50118250-CE	18	25cm	10/20mm	5
N50120060-CE	20	6cm	10/20mm	5
N50120100-CE	20	10cm	10/20mm	5
N50120130-CE	20	13cm	10/20mm	5
N50120160-CE	20	16cm	10/20mm	5
N50120200-CE	20	20cm	10/20mm	5
N50120250-CE	20	25cm	10/20mm	5

M·Biopsy® Semi-automatic Disposable Biopsy Needle Set

Ref. no.	Gauge	Length	Throw Length	Matching Coaxial Needle			Box
				Gauge	Length	With Blunt Needle	
NS5011410-CE	14	10cm	10/20mm	13	4,8cm	Yes	5
NS5011413-CE	14	13cm	10/20mm	13	7,8cm	Yes	5
NS5011416-CE	14	16cm	10/20mm	13	10,8cm	Yes	5
NS5011610-CE	16	10cm	10/20mm	15	4,8cm	Yes	5
NS5011613-CE	16	13cm	10/20mm	15	7,8cm	Yes	5
NS5011616-CE	16	16cm	10/20mm	15	10,8cm	Yes	5
NS5011620-CE	16	20cm	10/20mm	15	14,8cm	Yes	5
NS5011810-CE	18	10cm	10/20mm	17	4,8cm	Yes	5
NS5011813-CE	18	13cm	10/20mm	17	7,8cm	Yes	5
NS5011816-CE	18	16cm	10/20mm	17	10,8cm	Yes	5
NS5011820-CE	18	20cm	10/20mm	17	14,8cm	Yes	5
NS5012010-CE	20	10cm	10/20mm	19	4,8cm	Yes	5
NS5012013-CE	20	13cm	10/20mm	19	7,8cm	Yes	5
NS5012016-CE	20	16cm	10/20mm	19	10,8cm	Yes	5
NS5012020-CE	20	20cm	10/20mm	19	14,8cm	Yes	5
NS5012025-CE	20	25cm	10/20mm	19	17,8cm	Yes	5





mermaid medical®

NORDICS

Frydensbergvej 25 - DK-3660 Stenløse

customer@mermaidmedical.com

Phone: +45 47 10 85 70

Fax: +45 47 10 85 71

ITALY

customer@mermaidmedical.it

Phone: +39 010 868 16 13

NETHERLANDS

info@mermaidmedical.nl

Phone: +31 850 60 8160

SPAIN

customer@mermaidmedical.es

Phone: +34 925 717 378

UNITED KINGDOM

customer@mermaidmedical.co.uk

Phone: +44 1704 560493