



Ultra-Pro e™ Disposable Replacement Kits

The Ultra-Pro e™ is designed with simplicity in mind. It features an open-channel design with infinite angle capabilities, ensuring accurate in-plane positioning during procedures where access is difficult.

Features:

- Achieve faster, more accurate needle placement
- Accurate in-plane positioning with infinite angle options
- Compatible with all existing Ultra-Pro style brackets, assuring greater flexibility for puncture procedures.



The open channel design allows users an efficient needle release for viewing anatomy in both the short and long axis.



Ordering Information

Ultra-Pro e™ In-Plane Ultrasound Needle Guides

| Ref no. | Descriptions | Qty/box |
|----------|--|---------|
| 610-1082 | Sterile Ultra-Pro e™ 16G (guide only) | 24 |
| 610-1083 | Sterile Ultra-Pro e™ 16G with 14 x 91.5cm telescopically folded CIV-Flex™ cover | 24 |
| 610-1089 | Sterile Ultra-Pro e™ 14G with 14x91.5cm telescopically-folded CIV-Flex™ cover | 24 |
| 610-1133 | Sterile Ultra-Pro e™ 18G (guide only) | 24 |
| 610-1134 | Sterile Ultra-Pro e™ 18G with 14 x 91.5cm telescopically- folded CIV-Flex™ cover | 24 |
| 610-1137 | Sterile Ultra-Pro e™ 20G with 14 x 91.5cm telescopically- folded CIV-Flex™ cover | 24 |

MM 02-1000-06 ENG Rev. 02.

Contact Details, Nordics

customer@mermaidmedical.com
mermaidmedical.com

Contact Details, NL

info@mermaidmedical.nl
mermaidmedical.com

Contact Details, UK

customer@mermaidmedical.co.uk ·
mermaidmedical.com