




Stereotix-2 Needle Guides



Needle guides without slot for wire positioning. Designed for accurate placement of the needle during stereotaxi. All needle guides are compatible with the acquired software from GE Healthcare.

Features:

- Provides horizontal or lateral approach for core biopsy and localization
- Compatible with Senovision, Stereotix DS, Senographe 200D and DMR systems
- Packaged 2 guides per package, 5 packages per case; sold per case
- Sterile item



Accurate placement of needles with horizontal or lateral approach

Ordering Information

Stereotix-2 Needle Guides/Bushings

Ref no.	Description	Hole Diameter mm	Gauge	Qty/box
E63101DL	Stereotix-2 Needle Guide	0.9	21	10
E63101DM	Stereotix-2 Needle Guide	1	20	10
E63101DP	Stereotix-2 Needle Guide	1.4	18	10
E63101DS	Stereotix-2 Needle Guide	2.25	14	10
E63101EB	Stereotix-2 Needle Guide	1	20	10
E63101EC	Stereotix-2 Needle Guide	1.2	19	10

Contact Details, Nordics

customer@mermaidmedical.com
mermaidmedical.com

Contact Details, NL

info@mermaidmedical.nl
mermaidmedical.com

Contact Details, UK

customer@mermaidmedical.co.uk ·
mermaidmedical.com